



ADVANCED VOICE COMMANDS

- Say **“Repeat 1, 1, 2 three times”** to chart repeated numbers more quickly
- Enable Speak Back PD while learning or practicing with headphones to have Bola speak your last pocket depth command back to you to ensure accuracy
- Enable Speak Back All while learning or practicing with headphones to have Bola speak your last command back to you to ensure accuracy

SPOT CHARTING

Spot charting allows you to quickly mark a different tooth without moving your cursor position away from the current tooth you are charting. The following commands work with both Bleeding and Suppuration:

- Say **“Bleeding on distal lingual of tooth number five”** to spot chart for bleeding
- Say **“Suppuration on distal lingual of tooth number five”** to spot chart for suppuration

SEXTANTS, QUADRANTS, AND FULL MOUTH

For severe cases, you may need to apply bleeding or suppuration to sections of the mouth or the full mouth. The following commands work with both Bleeding and Suppuration:

- Say **“Bleeding on sextant three”** to mark a sextant
- Say **“Suppuration on quadrant one”** to mark a quadrant

FURCATION

Furcation is rare and only applies to certain teeth. Specifically teeth 1-3, 14-16, 17-19, and 30-32.

- Say **“Furcation grade two”** to mark a furcation on the current tooth
- Say **“Furcation grade two on tooth number 3”** to spot chart a furcation on a specific tooth

SUPPURATION

- Say **“Suppuration”** to mark the buccal of the current tooth
- Say **“Suppuration all on tooth number 5”** to mark all positions of a specific tooth
- Say **“Suppuration mesial and distal”** to mark specific positions on the current tooth

MOBILITY

Our mobility command allows you to mark a tooth (from 1-4) for mobility.

- Say **“Mobility plus two”** to mark the current tooth with a +2 measurement
- Say **“Mobility tooth number 5 plus three”** to mark a specific tooth with a +3 measurement

ADDITIONAL RESOURCES

[Perio Fundamentals \(Level 1\)](#)